



TZ Medical
Sparked by your ideas

STERILE DEFIB PADS




TZ Medical Sterile Defibrillation Electrodes are constructed to adhere and conform to the patient's body, providing superior current distribution, which results in better clinical outcomes.

Our multi-function sterile electrodes are ideal for a variety of procedures, including TAVR, AVR, Mini-Maze, ICD Implants, Laser Lead Extraction and much more. They are designed for defibrillation, synchronized cardioversion, external pacing and ECG monitoring. Compatible with biphasic and monophasic defibrillators, our sterile electrodes offer direct connection to Zoll® defibrillators. We also offer adapters, allowing multiple brands of defibrillators to use the same electrodes.



Radiotransparency



TZ

Competitor

Features

- Lead wire lengths available in 42", 60" and 120" for maximum convenience and utility
- Hydrogel maintains excellent adhesion to the patient's skin
- Unique edge effect technology disperses energy uniformly, reducing burns
- Radiotransparent and radiotranslucent options available

FIND OUT MORE



tzmedical.com



800.944.0187




info@tzmedical.com

13622 R2.3


Pad Placement Guide



TZ Medical
Sparked by your ideas



Lead Wire Length



Medtronic® M1 Connector



Philips® P2 Connector



Universal T1 Connector



Zoll® Z1 Connector



Product Description	Size	Medtronic® ¹	Philips® ²	Universal ³	Zoll® ⁴
Adult Radiotranslucent, sterile	A	P-121-T1 PAD 1111 Adapter*	P-121-T1 PAD 109 Adapter*	P-121-T1	P-121-T1 PAD 101PD Adapter*
Adult Radiotransparent pads and leadwires, large posterior pad, sterile	A/B	P-224-T1 PAD 1111 Adapter*	P-224-T1 PAD 109 Adapter*	P-224-T1	P-224-Z1
Pediatric Radiotranslucent pads, sterile	C	P-122-T1 PAD 1111 Adapter*	P-122-T1 PAD 109 Adapter*	P-122-T1	P-122-T1 PAD 101PD Adapter*

Pad Dimensions:

A=6"x 4.25" B=6"x9" C=4.25"x2.876"

For Pediatric Use: ≥10kg (22lbs)
AED or AED Mode: ≥25kg (55lbs) ≥8 years old

***Adapters are reusable**

- 1 – M1 Medtronic models are comparable to Medtronic® Physio-Control QUIK-COMBO™ electrodes.
- 2 – P2 Philips models are comparable to Philips® gray flat-style double barrel electrodes.
- 3 – T1 Anderson models allow for multiple brand of defibrillators to use the same electrodes. (adapter needed)
- 4 – Z1 Zoll Models are comparable to Zoll® multifunction electrodes.

FIND OUT MORE



tzmedical.com



800.944.0187



info@tzmedical.com

13622 R2.3