



TZ Medical
ClearBeat Monitoring Service



TRIDENT™ PRO

Holter • Extended Holter • Event • MCT



One Trusted Partner. All of Your Cardiac Monitoring Needs.

OUR DEVICE

INNOVATION MEETS SIMPLICITY

Trident™ Pro cardiac monitors combine decades of research and innovation to deliver a solution designed to streamline remote ECG recording from start to finish. Built-in telecom and rechargeable lithium battery deliver the simplicity of a single-piece device with options for remote or local configuration and data transmission.



Records up to 30 days of continuous ECG data



Large LCD color touchscreen with an on-board digital diary



Lithium battery may be charged while in use



Compact and discreet device with water resistant case

FEWER COMPONENTS. SMARTER MONITORING.

Trident™ Pro monitors are equipped with a large LED touch screen and on-board patient diary that provides a user-friendly experience for manually recorded patient events. Three channels of ECG data are acquired from surface electrodes adhered to the body, and 8G of internal memory with advanced detection algorithm offer a versatile solution to capture, store, and transmit every beat during Holter, Extended Holter, Event, and Mobile Cardiac Telemetry studies for up to 30 days.

MONITORS	DESCRIPTION
Trident™ 40	Digital Holter Monitor for Holter only studies
Trident™ 40L	Digital 4-in-1 monitor for Holter, Extended Holter, MCT, and Event studies, LTE wireless coverage

One Trusted Partner. All of Your Cardiac Monitoring Needs.

OUR SERVICE




Fast & Reliable Arrhythmia Detection

Our certified cardiac technicians provide accurate and reliable analysis of Holter, Extended Holter, Event, and Mobile Cardiac Telemetry data 24/7, rapid notification of emergent events, and industry leading report turnaround times.

Customer Support

Our friendly customer service representatives are available to serve you and your patients 24/7. No lengthy teleprompt, no unnecessary hold-times. Just live patient and provider support when you need it, every time.



- Simplified enrollments
- One-step device configuration & reset
- Industry leading report turnaround times
- Review & e-sign reports from any web-connected device or facility EMR system
- Seamless transition from Holter to Event or MCT

Contact us at info@clearbeat.com to learn more!

WHY CHOOSE US

ClearBeat is a comprehensive solution for your practice's Holter, Extended Holter, Cardiac Event, and Mobile Cardiac Telemetry needs. One trusted partner to advance your patients from symptom, to diagnosis, to treatment as quickly, simply, and cost-effectively as possible.

WHAT SETS US APART

Intuitive, durable, and state-of-the-art monitoring devices coupled with industry-leading report turnaround times, a simple and efficient workflow, and dedication to service after the sale make ClearBeat a clear-cut choice for your ambulatory cardiac monitoring needs.



TZ Medical Inc.
17750 SW Upper Boones Ferry Rd. #150
Portland, OR 97224
Office: 800.944.0187 Fax: 503.639.0239
www.tzmedical.com
18646 R1.2